µStat 8000 Multi Potentiostat/Galvanostat



()1 Ref. STAT8000





With a size of only 22x20x7 cm, **µStat 8000** is an equipment that includes **8 channels** that can act at the same time as **8 independent potentiostats/galvanostats**; it also includes **one multichannel** that can act as a potentiostat where up to 8 working electrodes share an auxiliary and a reference electrode.

With µStat 8000 users can perform up to 8 different electrochemical techniques at the same time; or carry out the study of one technique's parameter in just one step by applying the same electrochemical technique in several channels but selecting different values for the parameter under study. These are just examples of the enormous capabilities that this instrument offers.

μStat 8000 can be applied for Voltammetric, Amperometric or Potentiometric measurements, including 20 electroanalytical techniques.

The multi potentiostat/galvanostat is **Li-ion Battery powered** (DC charger adaptor also compatible), and can be easily connected to a PC via USB or **through Wireless connection**.

μStat 8000 is controlled by the included **software "DropView 8400"** which allows plotting of the measurements and performing the analysis of results. DropView software provides powerful functions such as experimental control, graphs or file handling, among others.

Available techniques:

<u>POTENTIOSTAT</u>

Voltammetry

LSV Linear Sweep Voltammetry CV Cyclic Voltammetry **SWV** Square Wave Voltammetry DPV Differential Pulse Voltammetry NPV Normal Pulse Voltammetry **NDPV** Differential Normal Pulse Voltammetry AC Voltammetry (only EC mode) ACV Linear Polarization Resistance LPR

<u>Amperometry</u>

AD Amperometric Detection

ZRA Zero Resistance Amperometry

FA Fast Amperometry (t_{int} < 0.1 s)

PAD Pulsed Amperometric Detection

COUL Coulometric Detection

<u>GALVANOSTAT</u>

LSP	Linear Sweep Potentiometry		
CP	Cyclic Potentiometry		
PD	Potentiometric Detection (galvanostatic)		
FP	Fast Potentiometry (t _{int} < 0.1s)		
ZCP	Zero Current Potentiometry		
PSAG	Potentiometric Stripping Analysis (galvanostatic)		
PSAF	Potentiometric Stripping Analysis (faradaic)		

Instrument Specifications					
Power	Li-ion Battery (6150 mAh) USB DC charger adaptor compatible (5 V, 15 W)				
PC interface	Wireless connection USB				
Operating modes	8x 1 Channel Potentiostat/Galvanostat 1x 8 Channel Potentiostat				
DC-Potential range	±4 V				
Current ranges (potentiostat)	± 1 nA to ± 100 mA (9 ranges)				
Maximum measurable current	±80 mA				
Potential ranges (galvanostat)	± 100 mV, ± 1 V (2 ranges)				
Applied Potential Resolution	1 mV				
Measured Current Resolution	0.025 % of current range (1 pA on lowest current range)				
Applied Current Resolution	0.1 % of current output range				
Measured Potential Resolution	0.012 % of potential range				
Potential Accuracy	±0.2 %				
Current Accuracy	≤0.5 % of current range at 100 nA to 1 mA ≤1 % of current range at 10 mA to 100 mA				
External inputs/outputs	 5 Digital Input/Output pins [PIO 1, PIO 2, PIO 3, PIO 4, PIO 5] 3 Analog Inputs multiplexing PIO 1, PIO 2, PIO 3 2 Analog Outputs (configurable I-out or E-out) 				
Indicators	LCD display in front panel				
Dimensions Weight	22.2 cm x 20.5 cm x 7.5 cm (L x W x H) 1.6 kg				

Control Specifications						
General Pretreatment	Conditioning stage duration: Deposition stage duration: Equilibration stage duration:	0 – 1300 s 0 – 1300 s 0 – 1300 s				
General Parameters	Begin, End, Base, Vertex potentials: Step potential: Pulse potential: Scan rate:	-4 V to +4 V 1 mV to 500 mV 1 mV to 250 mV 1 ms up to 1.3 s per step				
Specific Parameters	SWV	Frequency: Amplitude:	1 Hz to 400 Hz 1 mV to 250 mV			
	DPV, NPV, NDP	Modulation time: Pulse time:	1 ms to 1300 ms 1 ms to 1300 ms			
	ACV	Frequency: Amplitude:	2 Hz to 250 Hz 5 mV to 250 mV (RMS)			
	LPR	dE/dt lim: tmax OCP: tprecond:	-1 μV/s to 1000 μV/s 5 s to 6550 s 0 s to 1300 s			
	Chrono. Methods (AD, PD, ZCP, ZRA, COUL)	Interval time: Run time:	0.1 s to 1300 s Hours (65000 points)			
	Fast Chrono. Methods (FA, FP)	Interval time: Run time:	1 ms to 1300 ms Hours (65000 points)			
	PAD	Pulse time: Interval time: Run time:	1 ms to 1300 ms 10 ms to 1300 ms Hours (65000 points)			
	PSA	Potential limit:	±4 V			

Specifications are subject to change without previous notice

 $www.metrohm\hbox{-}dropsens.com$

