



DS APPLICATION AREA: FUNDAMENTALS

K_{cell} determination through impedance measurements

Summary

Conductivity of a solution (κ) can be related with a constant parameter of the sensor, defined as cell constant (K_{cell}), according to the equation:

$$\kappa = \frac{K_{cell}}{R}$$

Where R is the resistance measured.

In general, the higher the conductivity of the solution, the greater the K_{cell} required.

The cell constant value is related to electrode geometrical parameters (such as digit length,

number of digit pairs, gap between digits...). This constant can be determined experimentally¹. The cell constant is defined as:

$$K_{cell} = \frac{R_{sol}}{\rho_{sol}}$$

Where R_{sol} is the resistance of the medium and ρ_{sol} is electrolyte resistivity.

The cell constant could be determined by immersing the electrode in a standard solution, generally potassium chloride, whose conductivity at different

temperatures and different concentrations is known. Impedance measurements allow to determine R_{sol} at frequencies where phase angle is 0 or close to 0 (Bode plots).

Apparatus and accessories:

Impedance measurements are carried out with an Autolab PGSTAT 204 controlled by NOVA 1.10 software. Impedance spectra were taken in the range from 100 Hz to 1×10^6 Hz, RMS 10mV (peak-to-peak amplitude). There is an equilibration time before each measurement of 5 s. Measurements were performed with IDEs, multi-interdigitated electrodes and SPEs of different geometry and materials.

K_{cell} is calculated as the slope from the plot Impedance ($|Z|$) vs Solution resistance (ρ).

Reagents:

KCl solutions are made from dilutions of a conductivity standard 1.0M (111,8 mS/cm). The conductivity from KCl solutions is calculated by Kohlrausch equation ($EC = \sum(c_i \times f_i)$; $c_i = \text{mg/L}$ and f_i is the conductivity factor).

Method:

The impedance spectra of each sensor immersed in the different KCl solutions were registered. All measurements were done at 25-26°C. $|Z|$ (or solution resistance) is plotted versus the inverse of the conductivity (or resistivity) at a frequency where ϕ or angle phase is close to 0 or has minimum value. The slope of these curves is the electrodes cell constant.

Results:

Electrode REF.	Frequency (KHz)	Conductivity range ($\mu\text{S cm}^{-1}$)	K_{cell} (cm^{-1})
G-IDEAU5	200	84 - 12880	0.006
G-IDEAG5	200	84 - 1413	0.006
G-IDEPT5	100	84 - 12880	0.006
G-IDEAU5	500	84 - 12880	0.006
G-IDEAU10	200	84 - 12880	0.012
G-IDEPT10	200	84 - 12880	0.012
G-IDECONAU10	200	84 - 12880	0.013
G-IDECONPT10	200	84 - 1413	0.013
P-IDEAG50	200	84 - 12800	0.059
P-IDEITO50	200	84 - 1413	0.059
PW-IDEPD50	200	84 - 12800	0.059

PW-IDEAU50	32	84 - 12800	0.059
PW-IDEPD100	200	84 - 12800	0.124
P-IDEITO100	25	84 - 1413	0.124
P-IDEAG100	25	84 - 12800	0.124
PW-IDEAU100	100	84 - 12800	0.124
P-IDEAU100	100	84 - 12800	0.124
IDEAU200-HPT-WB	200	84 - 12800	0.220
IDEAU200	200	84 - 12800	0.220
G-IDE222	400	84 - 12800	0.107
G-IDE555	400	84 - 12800	0.107
PW-4XIDEAU20	200	84 - 12800	0.300
PW-4XIDEAU30	200	84 - 12800	0.467
PW-4XIDEAU40	200	84 - 12800	0.600
PW-4XIDEAU50	200	84 - 12800	0.700
PW-4XIDEAU60	200	84 - 12800	1.050
PW-4XIDEAU70	200	84 - 12800	1.050
TLFCL222AT	50119	84 - 12800	10
X1110	50119	84 - 12800	2.710

References:

¹ De la Rica, R.; Fernández-Sánchez, C.; Baldi, A. *Electrochem. Commun.* **2006**, 8, 1239.