

# **DROPSENS**

## *Flow-Cell for Screen-Printed Electrodes*

(ref. FLWCL)



**DropSens Flow-Cell (ref. FLWCL) is a methacrylate wall-jet flow-cell for FIA, with an innovative open-close system (no screws needed) and an easy sensors replacement**



**An o-ring limits the volume of the electrochemical cell**



**It can be used with all DropSens screen-printed electrodes**





(+34) 985 27 76 85

Edificio Severo Ochoa, Local 5 - El Cristo  
33006 Oviedo (Asturias). Spain


[www.dropsens.com](http://www.dropsens.com)

This cell is designed to obtain an inlet flow perpendicular to the electrode's surface



Due to the transparency of the material, the possible presence of bubbles is easily detectable

The use of Cable Connector (ref. CAC) allows an easy connection to any potentiostat



**DropSens, S.L.**

Edif. Severo Ochoa, Local 5 - El Cristo  
33006 Oviedo (Asturias) - Spain

Phones: (+34) 985277685 / 653525278

FAX: (+34) 985277685

[info@dropsens.com](mailto:info@dropsens.com)

[www.dropsens.com](http://www.dropsens.com)